

Delikit

KOREA NO.1

VPS Dental Impression Materials



Delikit Aqua-Kicker Light Body

Technical data

Working time	1' 30"
Setting time (in the mouth)	4' 00"

Advantages

1. Impression of inlay, onlay, crown, and bridge preparation
2. Excellent operability, compare with competitor



50ml

Delikit Impre-Kicker Light Body

Technical data

Working time	1' 30"
Setting time (in the mouth)	4' 00"

Advantages

1. Impression of inlay, onlay, crown, and bridge preparation
2. Excellent operability, compare with competitor



50ml

Delikit Light Body - Regular

Technical data

Working time	1' 30"
Setting time (in the mouth)	4' 00"

Advantages

1. Nice flowability by air blowing
2. Thixotropic
3. Excellent wetting of humid areas
4. High tear strength
5. Working time is secured sufficiently by Snap-set technology.
6. High dimensional stability
7. Maximum precision in detail reproduction



50ml

Delikit Light Body - Extra

Technical data

Working time	1' 30"
Setting time (in the mouth)	4' 00"

Advantages

1. Suitable for inlay, onlay, single crown
2. High tear strength
3. Excellent elastic recovery
4. Possibility of casting the model immediately after disinfection
5. Maximum precision in detail reproduction
6. High dimensional stability



50ml

Delikit Light Body - Fast

Technical data

Working time	1' 15"
Setting time (in the mouth)	2' 30"


Advantages

1. Suitable for small case ex) inlay, onlay, single crown
2. Quick set
3. High tear strength
4. Excellent reproduction of details
5. Very good mechanical properties
6. Excellent elastic recovery
7. High dimensional stability
8. No release of toxic substances during and after curing



50ml

Secret of error-free




**combination of accessibility,
controllability,
speed,
and individuality.**

INTRODUCE

■ Mirco Light Body: Single use syringe!

It is able to be accessed to prepared tooth as close as possible. Thanks to long and thin tip, it is inserted to proximal area without bubbles. Pencil like Syringe, and it is able to reduce working time. Therefore, errors due to initial setting in the mouth are avoided. Errors due to overall defextiveness of cartidge are removes from individual syringe.





Secret i : Light body Fast

- 1.7ml (for single case)
- W.T : 0'40" S.T : 2'00"



Secret ii : Light body

- 2.1ml (for single case)
- W.T : 1'30" S.T : 3'30"




Secret iii : Light body

- 2.7ml (for bridge)
- W.T : 1'30" S.T : 3'30"

HOW TO USE



Step1. Remove the tab



Step2. Push it in



Step3. Ready to use

Delikit Heavy Body - *Regular*

Technical data

Working time	1' 30"
Setting time (in the mouth)	4' 00"

Advantages

1. Easy to remove tray from the mouth
2. Choose your preferred technique-laminar
3. Thixotropic
4. No transformation by Dimensional stability
5. Excellent elastic recovery
6. High dimensional stability



50ml

Delikit Heavy Body - *Fast*

Technical data

Working time	1' 15"
Setting time (in the mouth)	2' 30"

Advantages

1. Inlay consistency for one step technique
2. Quick set
3. Suitable for small case
ex) inlay, onlay, single crown
4. good for disposable bite-tray



50ml

Delikit Putty - *Regular*

Technical data

Working time	1' 30"
Setting time (in the mouth)	3' 00"

Advantages

1. Short setting time in the mouth
2. Non-greasy, satiny smooth, Extremely easy to mix
3. No bubbles
4. Very good mechanical properties
5. Excellent reproduction of details
6. Hydrophilic



600ml

Delikit Putty - *Fast*

Technical data

Working time	1' 30"
Setting time (in the mouth)	2' 30"

Advantages

1. Suitable for inlay, onlay, single crown
2. Quick set



600ml

Bite registration material

Technical data

Working time	20"
Setting time (in the mouth)	50"

Advantages

1. Easy handling and bubble-free
2. Positioning control
3. It won't compress or flex
4. Thixotropic for easy application
5. Flowable: it avoids pressure and slipping during positioning
6. Bubble-free
7. Easy trim



50ml

Monophase/ Heavy body



380ml

Technical data

Working time	1' 30"
Setting time (in the mouth)	3' 30"

Advantages

1. Monophase material to be used both on tray and around the preparation for excellent direct reproduction
2. Ideal consistency for single-stage
3. Excellent elastic recovery
4. High dimensional stability
5. Possibility of casting the model immediately after disinfection

Dispenser 1:1/2:1

to be used with 50ml cartridges



Bite-tray



Posterior

Sideless
Posterior

Quadrant

Full arch

Anterior

Packaging

Posterior	50 ea/Pack
Sideless posterior	50 ea/Pack
Quadrant	30 ea/Pack
Full arch	25 ea/Pack
Anterior	30 ea/Pack

Mixing tip



**Mixing tip
(Light body)**



**Mixing tip
(Heavy body)**



**Mixing tip
(Jumbo)**



**Mixing tip
(Secret)
2.1/2.7ml**

Packaging

Mixing tips (Light body)	48 tips with 20 intraoral tip/Pack
Mixing tips (Heavy body)	48 tips/Pack
Mixing tips (Jumbo)	50 tips/Pack
Mixing tips (Secret) 2.1/2.7ml	50 tips/Pack

Overseas Certification



CE Certification



ISO 13485 Certification



FDA Certification

Hong Kong
Exclusive Distributor
Advance Dental Consulting Ltd.
www.adchk.com
info@adchk.com
Tel: 2384 7666 Fax: 2388 0960

